



HRM™

high-resolution melt



HRM APPLICATIONS

Andrea A Tesoriero



www.corbettlifescience.com

High Resolution Melt (HRM)

detailed monitoring of changes in fluorescence as a PCR products are melted

- ⊕ PCR amplify an amplicon with two primers and an intercalation dye and then melt (dissociate) the product – double stranded to single stranded
- ⊕ Can detect even a single base pair change

HRM Specifications

⊕ **Instrument requires:**

- ◆ high-intensity + high sensitivity optics
- ◆ high-speed data capture
- ◆ very precise temperature control and resolution

Rotor-Gene™ 6000


The world's only real-time rotary thermo-optical analyzer

An ideal platform for HRM



Rotor-Gene™ 6000


Cross-section of rotary optics



Thermal Control ± 0.01


Rotor-Gene™ 6000

Heating mechanism



Rotor-Gene™ 6000

Cooling mechanism




HRM Profile

Quick Start

1. Rotor Selection | 2. Confirm Profile

New Open Save As Help

The run will take approximately 142 minute(s) to complete. The graph below represents the run to be performed :



Click on a cycle below to modify it :

Hold	Insert after...
Cycling	Insert before...
Hi-Res Melt	Remove


Ramp from 80 degrees to 90 degrees,
Rising by 0.02 degree(s) each step,
Wait for 90 seconds of pre-melt conditioning on first step,
Wait for 2 seconds for each step afterwards.
Acquire to Hi-Res Melt A on HRM

Gain Optimisation
 Optimise gain before melt on all tubes.
The gain giving the highest fluorescence less than 70 will be selected.

< Back Save Template **Start Run** Cancel


0.02°C

Data Acquisition




- ⊕ Melt curves normalized—by selecting linear regions before & after the melt transition
- ⊕ Two regions defined—upper 100% double stranded / lower single stranded baseline

Normalized Data: Homoduplexes C or T




Homozygotes represented by a single base change are differentiated by a difference in T_m

Heteroduplex C>T



Heterozygotes form heteroduplexes, which result in lower melt curves and two transition points

Difference Graphs



- ⊕ Difference graph displays each sample plot subtracted from a given genotype control plot
- ⊕ Calculates a percentage confidence relative to a known genotype

HRM Summary

- ⊕ Simple, fast, cost effective method for gene scanning and detecting a single-base change in your sample
- ⊕ Closed system
- ⊕ NO labeled probes
- ⊕ Has excellent sensitivity and specificity is capable of detecting BOTH heterozygous and homozygous changes
- ⊕ Costs less than competing technologies
- ⊕ Sequence directly off the product—sample not consumed
- ⊕ Auto call software
- ⊕ Scanning & genotyping can be done simultaneously in the same reaction

SNP's


Detection of alpha-thalassemia-1 Southeast Asian type using real-time gap-PCR with SYBR Green 1 and high Resolution Melting analysis. Pornprasert S, Phusua A, Suanta S, Saetung R, Sanguansermsri T, *Eur J of Haematol.* 2008.

Direct Genotyping of Single Nucleotide Polymorphisms in Methyl Metabolism Genes Using Probe-Free High-Resolution Melting Analysis. L S. Kristensen and A. Dobrovic *Cancer Epidemiol Biomarkers Prev* 2008;17(5) May 2008

High-resolution melt analysis for the detection of a mutation associated with permethrin resistance in a population of scabies mites C Pasay, L Arlian, M. Morgan, D Vyszynski-Moher, A. Rose, D Holt, S. Walton and JMcCarthy *Medical and Veterinary Entomology* (2008) **22, 82–88**

Genotyping Class 4 SNP

SNP Class	Base Change	Typical Tm Shift	Rarity (in humans)
1	C/T and G/A	Large >0.5oC ↓ Very Small >0.2oC	64%
2	C/A and G/T		20%
3	C/G		9%
4	A/T		7%





SNP classes as described by Venter *et al* 2002

Example of a class 4 SNP on the Rotor Gene (MCT A1470T)
The rarest and most difficult SNP to discriminate.


Class 4 SNP using Sybr Green 1

- SYBR can be used at non-saturating concentrations, for the most demanding HRM applications.



- Extra sensitivity of the Rotor-Gene may be a major reason we do not require saturating dye levels


Adjacent SNP's Homozygosis'



Homozygotes represented by a 2bp base changes are differentiated by a difference in T_m by about 0.3deg


Data Courtesy of Howard Florey Institute, J.Rubio, C. Jenson

Adjacent SNP's Heterozygous



Heterozygotes form heteroduplexes, which result in lower melt curves and multiple transition points

Nomalised Data



Sequences were performed directly from HRM product, variants are 10bp apart

Complex Melts -HLA Compatability

Relationship	04	11
Patient	●	●
Maternal Aunt 1	●	●
Paternal Aunt		●
Maternal Aunt 2		●


2 days

20 probes

SS sequencing

04	11
0403	1143
0401	1104
	X
	X

2 days




No match

Red Cross-Melbourne

Mutation Detection

Genomic Mutations in LDLR gene- 23 fragments



Atherosclerosis 194 (2007) 279–286. **Development of an affordable, sensitive and rapid screening method for mutation detection in UK FH subjects** C Hubbart, R Whittall, M Scartezini, S Humphries

White *et al.* 2006 report



Evaluation of High Resolution Melt Analysis for Mutation Scanning


Helen White, GR Taylor, GL Potts, NCP Cross, C Taylor

National Genetics Reference Lab (Wessex), Salisbury
Regional Genetics Service and CR-UK Mutation
Detection Facility, Leeds

<http://www.ngrl.org.uk/Wessex/downloads.htm>

White *et al.* 2006 report


Amplicons and Mutations Analysed I



White *et al.* 2006 report

Amplicons and Mutations Analysed II


hMLH1 Exon 1 (193bp, 57% GC Rich)




hMLH1 Exon 7 (139bp, 37% GC Rich)



hMLH1 Exon 13 (277bp, 44% GC Rich)



hMSH2 Exon 10 (249bp, 34% GC Rich)



White *et al.* 2006 report

Summary


- 624 samples analysed in total
- 11 amplicons tested: 139bp – 449bp with GC contents of 22-79%
- Mutations included all possible point mutation base substitutions and 1 and 2bp insertions and deletions
- The same PCR reaction was analysed using the HR-1 and 384 well LightScanner (Idaho Technology) and RotorGene 6000 (Corbett Research)

Sensitivity and Specificity of mutation detection for each platform were:


	Sensitivity	Specificity
RotorGene 6000	100.0	95.3
HR-1	98.4	95.0
LightScanner 384 well (High)	99.0	88.0
LightScanner 384 well (Normal)	83.9	95.3

DNA Quality—late amplification


White *et al.* 2006 report



DNA quality-Multiplex 100, 200, 300, 400 and 600 bp product




Amplification of 193 bp product



Sensitivity


“Spiking” Experiments

-applicable for acquired mutations and pooling samples




“Spiking” Experiments

-applicable for acquired mutations and pooling samples




“Spiking” Experiments

-applicable for acquired mutations and pooling samples



“Spiking” Experiments

-applicable for acquired mutations and pooling samples



Courtesy Royal Melbourne Hospital-Hematology

Acquired Mutation Detection

High resolution melting for mutation scanning of *TP53* exons 5-8

M Krypuy, A Ahmed, D Etemadmoghadam, S Hyland, Australian Ovarian Cancer Study Group, J Brenton, S Fox, A deFazio, D Bowtell and A Dobrovic *BMC Cancer* 2007, 7:168


High resolution melting analysis for the rapid and sensitive detection of mutations in clinical samples: *KRAS* codon 12 and 13 mutations in non-small cell lung cancer

M.Krypuy, G.Newnham, D.Thomas, M.Conron, A.Dobrovic *BMC Cancer*. 2006 Dec 21;6(1):295

Detection of the transforming *AKT1* mutation E17K in non-small cell lung cancer by high resolution melting . H Do, B.Solomon, P.Mitchell, S.Fox ,A.Dobrovic *BMC Research Notes* 2008, 1:14

Sensitivity

Krypuy *et al.* 2006 report




38 G>A
 Heterozygous
 92bp 5% cellular ratio
 1:19 dilution
 Represents 2.5% allele ratio
 1 out of 39 alleles

34 G>A
 Homozygous

35 G>C
 hypotriploid
 possible 2:1 mutant

Somatic Mutation Discovery KRAS


Exon 2 189 bp product 37% GC content-



- Detected small quantities of mutant DNA in a background of wildtype DNA
- Confirmed 9/30 positive clinical samples
- No homozygous spiking necessary
- The heteroduplex shape is not used to determine differences, rather the deviation from wildtype, as the shape is compromised due to the proportion of normal DNA

Difference Graph

Krypuy *et al.* 2006 report



This paper concludes that HRM will allow

*“confident screening for mutations in samples
that have at least 10% tumor cells”*

Krypuy *et al.* 2006 report

Species Identification

Classification of *Mycoplasma synoviae* strains using single-strand conformation polymorphism and high-resolution melting-curve analysis of the *vlhA* gene single-copy region Nathan Jeffery, Robin B. Gasser, Penelope A. Steer and Amir H. Noormoham *Appl Environ Microbiol.* 2007; 73: 3431-3436

High-Resolution DNA curve analysis of the clustered regularly interspaced short palindromic-repeat locus of *Campylobacter jejuni*. Price EP, Smith H, Huygens F, Giffard PM *Microbiology (2007), 153, 2679–2688*

Strain Sequence Alignment

```

Consensus      1          11          21          31          41          51          61          71          81          91
TACTATTAGCAGCTAGTCCAGTCCGCCATTGCCCTGCTGTTATAGCAATTTTCATGTGGTGATCAAACTCCAGCACCTGCTCCA-----ACACC
94027/10a     .....A.....
T2/3X        .....T.....
94041/12a     .....T.....
93148/23-22b .....T.....
K1968        .....A.....
MVU-1853     .....AA.....
YA           .....AA.....
F10-2AS     .....AA.....
MS-H        .....
94046/M1B-17a.....
94011/V-18d.....
94042/6a     .....
94029/1a     .....
K1723        .....
K1858        .....G.....
K1938        .....G.....A.....GCACCTACTCCA.....

Consensus     101         111         121         131         141         151         161         171         181         191
TGGAARCCCAATACTGATATCTCTCAAAAACCCAAATCCAGGAAA CCAGS ACT-----GA-----TAAT
94027/10a     .....C.....G.A.GATAATCTCTCAAAAACCCAAATCCAGGAAAACCCAG.....
T2/3X        .....C.....G.A.GATAATCTCTCAAAAACCCAAATCCAGGAAAACCCAG.....
94041/12a     .....C.....A.....GATAATCTCTCAAAAACCCAAATCCAGGAAAACCCAG.....
93148/23-22b .....C.....A.....GATAATCTCTCAAAAACCCAAATCCAGGAAAACCCAG.....
K1968        .....C.....A.....GATAATCTCTCAAAAACCCAAATCCAGGAAAACCCAG.....ACTGA.....
MVU-1853     .....T.....T.....CCAGGAAAATCCAGSSTACT.....
YA           .....T.....T.....CCAGGAAAATCCAGSSTACT.....
F10-2AS     .....
MS-H        .....
94046/M1B-17a.....
94011/V-18d.....
94042/6a     .....T.....
94029/1a     .....T.....
K1723        .....T.....
K1858        .....
K1938        .....

Consensus     201         211         221         231         241         251         261         271         281         291
CTCAARACCCAAATCCAGGAAA CCAGS GGINGGTACAGTTGAACCTGTAGAGGCTCTCTAAAACAGAGACTAAAACCCCTATTGATGCTTCAGCAGAAT
94027/10a     .....T.....G.A.....
T2/3X        .....T.....G.A.....
94041/12a     .....T.....G.A.....
93148/23-22b .....T.....G.A.....
K1968        C.....T.....A.....
MVU-1853     C.....C.....A.....
YA           C.....C.....A.....
F10-2AS     .....C.....A.....
MS-H        T.....T.....G.....
94046/M1B-17a T.....T.....G.....
94011/V-18d T.....T.....G.....
94042/6a     T.....T.....G.....
94029/1a     T.....T.....G.....
K1723        C.....T.....A.....
K1858        .....G.....
K1938        .....G.....A.....T.....A.....



Consensus     301         311         321         331         341         351         361         371         381         391
TATCAGATTTCAGTTAAAGGAGCATTAAAAGACACAGTTGANGCAACTACACACGAGCTCCAGCCGAGATTTAAAACCTAAGCAGANGCTCTTGTTC
94027/10a     .....C.....A.....
T2/3X        .....C.....A.....
94041/12a     .....C.....A.....
93148/23-22b .....C.....A.....
K1968        .....C.....A.....
MVU-1853     .....A.....
YA           .....A.....
F10-2AS     .....A.....
MS-H        .....
94046/M1B-17a .....
94011/V-18d .....
94042/6a     .....
94029/1a     .....
K1723        .....
K1858        .....
K1938        .....T.....

Consensus     401         411         421


```


SSCP

400bp PCR Product



Normalised HRM



~400 bp fragment, 35 different *M. synoviae* isolates

Epigenetics

Methylation-sensitive high resolution melting (MS-HRM):
A new approach for sensitive and
high-throughput assessment of methylation.

Tomasz K. Wojdacz, Alexander Dobrovic

Nucleic Acids Research Feb 2007

***“The simplicity and high reproducibility...
...makes MS-HRM the method of choice”***

What is methylation?

- ⊕ Addition of a methyl group to cytosine DNA in CpG island
- ⊕ Bisulphate treatment converts non methylated cytosine bases to uracil and methylated remain unaffected
- ⊕ Uracil converted to thymine by PCR




- ⊕ Methylation acts to “turn off” transcription
- ⊕ Markers for; Early events in cancer
Monitoring progression
MGMT predictive markers treatments

MGMT MS-HRM


Wojdacz *et al.* 2007 report

Effect of annealing temperature on sensitivity




60°C

Discrimination at higher end



62°C



63°C

Discrimination at lower end




increasing temperature


Low level methylation 0.1% MGMT sensitivity MS-HRM

Wojdacz *et al.* 2007 report

12 CpG's/ 109bp




5CpG's/ 94bp



Heterogeneous methylation

Wojdacz *et al.* 2007 report



MS-HRM is a sensitive, reliable,
fast and cost effective

“...*method of choice*”

Wojdacz *et al.* 2007 report

Genomic Imprinting

Rapid detection of methylation change at H19 in human imprinting disorders using methylation-sensitive high-resolution melting.

Wojdacz TK, Dobrovic A, Algar EM. *Hum Mutat.* 2008 May 12

Methylation-Sensitive High-Resolution Melting-Curve Analysis of the SNRPN Gene as a Diagnostic Screen for Prader-Willi and Angelman Syndromes. H White, V Hall, and N Cross. *Clinical chemistry* September 2007

Promoter Methylation

BRCA1 promoter methylation in peripheral blood DNA of mutation negative familial breast cancer patients with a BRCA1 tumour phenotype. C Snell, M Krypuy, E Wong, M Loughrey, A Dobrovic
Breast Cancer Research 2008

Applications

- ⊕ Detect small quantities of mutant DNA in background of wildtype DNA species
- ⊕ Important in somatically acquired mutations
- ⊕ Pooling samples—up to 10 samples
- ⊕ Simple for diseases that cause no heterogeneity-like Factor V Leiden, haemochromotosis, sickle cell anemia, etc
- ⊕ Appropriate for disorders where the mutational spectrum is wide—mismatch repair genes
- ⊕ Newly identified genes—little information
- ⊕ Species identification
- ⊕ Epigenetics




www.corbettlifescience.com



www.corbettlifescience.com

Singleplex HRM



Multiplex HRM

